

Microstructural Characterization of GaN LEDs by TEM


GaN LEDs are fabricated on substrates such as sapphire, but lattice mismatch causes dislocations to form. Using TEM, the types and densities of these dislocations were evaluated. Moreover, the MQW structure at the center of the bright region was analyzed using TEM and EDX.

Dislocation Analysis in GaN


The dislocation type is determined based on the intensity observed in TEM images under the two-beam diffraction condition.

The dislocation density is calculated using Ham's method. (The table on the right shows the results.)

Dislocation type	Dislocation density (1/cm ²)
Edge Dislocation (red)	6.8E+8
Screw Dislocation (yellow)	4.9E+7
Mixed Dislocation (blue)	6.6E+8
Total	1.4E+9



Lattice-image observation of the MQW structure



Connecting technology, information, and people to give shape to "willingness"

Contact us:
<https://iontc.co.jp/contact>

